INSTRUCTIONS FOR USE

UroDapter® Urological syringe adapter

Thank you for purchasing **UroDapter®**, which enables catheter-free instillation of the bladder.

Compared to the traditional way of bladder instillation – which involves a catheter – using **UroDapter®** has several advantages. The process is minimally invasive and entirely painless. A catheter, even if it is used in the most cautious way, can cause microlesions and be the source of catheter-associated infections. These issues are not experienced if **UroDapter®** is applied. Moreover, **UroDapter®** enables the treatment of the bladder and the urethra at the same time, which is challenging or impossible with a catheter.

INTENDED USE

Performing bladder instillation with any solution prescribed by a therapist without a catheter, in a safe, painless, and minimally invasive way.

Generally used solutions include substances for regenerating the inner mucosal surface of the bladder - the GAG-layer - (e.g., hyaluronic acid, chondroitin sulphate, heparin, DMSO, and any combination of them), local anesthetics (e.g. lidocaine); antibiotics; antiseptics; certain hormones, chemotherapeutic agents (mitomycin C, epirubicin, gemcitabine, valrubicin, docetaxel) or immunotherapeutic agents (BCG - Bacillus Calmette-Guérin) to treat bladder cancer, and agents to treat hemorrhagic cystitis.

INDICATIONS

Indications for which **UroDapter®** can be successfully applied are Interstitial Cystitis/Bladder Pain Syndrome (IC/BPS), recurrent Urinary Tract Infections (rUTIs), certain forms of Hemorrhagic Cystitis (Chemotherapy Induced Cystitis, Radiation Cystitis), Severe Overactive Bladder Syndrome (OAB), and local chemo- or immunotherapy of Bladder Cancer. Product could be used also with physiological saline, phage solutions, anti-inflammatory agents (NSAIDS, corticosteroids), beta-3 adrenergic agonists (mirabegron), anticholinergic agents (oxybutynin, tolterodine, solifenacin, trospium, darifenacin), and contrast agents for imaging. **UroDapter®** can be applied to male or female patients.

CONTRAINDICATIONS

Use of product is not advised on patients with known or suspected hypersensitivity to any of the materials employed.

INTENDED USERS

UroDapter® can be used by medical professionals who are expert in the field of bladder instillation and experienced in the use of the device.


UroDapter® can be used by any non-professional person providing they have been trained on how to use the device by a professional priorly.

MATERIALS USED

The product is made of medical grade thermoplastic elastomer (TPE).

INSTRUCTIONS FOR USE

1. Perform the instillation only in places with suitable hygiene level – poor hygienic environment may result in complications. Before starting the treatment, disinfect not only the urethral orifice but also the 3 cm surrounding area and the inner surface of the labia with any mucosal antiseptic solution.
2. Discard the product if the packaging is damaged. Unpack **UroDapter®** and attach its connective tail (4) to the syringe. After taking the product out of the packaging, inspect the product visually for mechanical damage – do not use the product if there is any present. After unpacking, use the device immediately. **UroDapter®** tail (4) can be attached to both Luer Slip and Luer Lock syringes easily. Although **UroDapter®** is sterile, it might be worth applying disinfectant agent on it, from the tip (1) to the sealing collar (2). This process can prevent the drift of the bacteria into the bladder, and inactivate the bacteria located in the first couple of millimeters of the urethra.

Parts of UroDapter®

1. **Radiused tip** – this enters the urethra (appr. 6–8mm deep).
 2. **Isolating collar** – this prevents leakage during the instillation.
 3. **Handling grip** – by this part **UroDapter®** can be held firmly.
 4. **Connecting tail** – this part is to be connected to any Luer Lock or Luer Slip syringe.
3. Always make it sure that **UroDapter®** is inserted into the urethra, and not into the vagina in case of female patients.
 4. For exposing the orifice of female patients, the fingers should lie as close to the orifice as they can be. The exposition has to happen by pulling the labia apart laterally, too (and not only upwards).
 5. The optimal instillation can be performed if the axes of the urethra and the syringe are aligned in parallel. In this case, the resistance can be experienced by using the syringe decreases characteristically.

6. Only the tip of **UroDapter®** should penetrate the urethra. The isolating collar (2) has to cover the orifice properly and be pushed against the nearby structures gently.
7. Start performing the instillation. Given that the isolating collar (2) is covering the orifice properly, the instilled solution raises the intra-urethral pressure, which opens the sphincter, thus the solution enters the bladder.
8. If the sphincter becomes constricted due to pain or stress, tell the patient to relax. Wait until the constrictions end: the speed of the instillation should be decreased or paused.
9. In case of leakage, a part of the solution does not reach the urethra, it flows away instead. The two most frequent situations in which leakage happens are applying too low or too high pressure with **UroDapter®**. Adjust the pressure accordingly; the ideal amount can be experienced as the sudden lowering in resistance of the plunger of the syringe.
10. In most cases, the leaking solution causes neither infections nor complications by entering the vagina or contacting the nearby structures. In case of doubt, consult the Instructions for Use or general description of the solution being instilled.

PRECAUTIONS AND WARNINGS

- **UroDapter®** is sterile; it is for single use only.
- Only a professional person can prescribe the use of **UroDapter®** for self-instillation. Use is allowed to professionals or priorly trained individuals.
- Before starting a self-therapy with **UroDapter®** always consult with your therapist about
 - the way how the self-treatment should be performed
 - the solution you have to instill into your bladder
 - the frequency and the total number of self-treatments.
- **UroDapter®** is designed for administering solutions into the bladder, through the orifice of the urethra. Never use **UroDapter®** for delivering any liquid into any other organ.
- **UroDapter®** is intended for instilling solution into the bladder. Do not use it for draining of urine or any other liquid.
- Do not reuse or resterilize **UroDapter®**.

DANGER!

- In case of reuse, risk of transmitting infections and occurrence of inflammation is present.
- Do not use the product if its packaging is damaged.
- After use, product has to be discarded in accordance with medical practice and local regulations.
- Keep **UroDapter®** away from children younger than 3 years, since it is a suffocation hazard for them.
- In case of performing an instillation with the help of a non-professional person, ensure that they have read

this Instructions for Use before starting the treatment, and have understood the instructions of your therapist fully.

- Ensure that you will not have to change posture during the self-treatment against your will; do not commit public indecency.
- In case of any side effects caused by the treatment or the instilled solution, contact your therapist immediately.
- In case of adverse effects caused directly by **UroDapter®**, inform Pharmico Medical Kft. on the incident.
- In case the product is not used in accordance with the instructions provided by Pharmico Medical Kft., the manufacturer or the distributor do not take any responsibility for any injuries or other damage that may arise.
- For any complications or adverse effects caused by the solution instilled with **UroDapter®** neither the manufacturer, nor the distributor take any responsibility.

STORAGE

If the product or its sterile packaging is damaged, the risk of transferring infections may be increased/introduced. Moisture or inappropriate temperature in the storage room can disrupt product sterility, thus, affect its safety and usability.

The product should be stored in its original packaging in a dry place, in a temperature range from 5 to 35°C (from 41 to 95°F).

The product should be protected from extensive heat or direct sunlight.

The sterility for a properly stored product is guaranteed for 5 years unless the packaging is opened or damaged.

For further information please visit our websites: www.urosistem.com, www.urodapter.com and <https://app.urodapter.com> where you can find a detailed Tips & Tricks guide of using **UroDapter®** for therapists.

This Instructions for Use is a subject of updates. The latest version of the document is always available on the www.urosistem.com webpage.

The product is patent protected.

MANUFACTURER:

Pharmico Medical Kft.
4032 Debrecen, Füredi út 98. Hungary
E-mail: boss@pharmicomedical.hu

DISTRIBUTOR AND PATENT OWNER:

Urosystem Zrt.
1137 Budapest Szent István park 26. Fsz.2. Hungary
E-mail: mail@urosistem.com

GRAPHIC SYMBOLS EXPLANATION

	Sterilized using irradiation
	Catalogue number
	Unique device identification
	Batch code
	Medical device
	Manufacturer
	Date of manufacturing
	Distributor
	Use by
	Single use
	Non-pyrogenic
	Keep dry
	Temperature limit
	Do not use if packaging is damaged
	PVC free
	Caution
	Consult instructions for use
	Latex free

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